

Phoenix 5000

NIR Feed Analyzer



Analyze feed ingredients, finished feed and by-products

9 30 Results in 30 seconds for moisture, fat, protein, fiber, starch and more

- Results in 30 seconds or less
- Monitor incoming ingredients
- > Determine least cost formulations
- Control moisture in driers and extruders
- > Verify finished product
- Transfer existing databases from other NIR instruments
- > User friendly software design
- Flexible data management

Next Generation NIR Analyzer

The Phoenix 5000 NIR Feed Analyzer is an ideal solution for feed producers that need highly accurate rapid results before, during and after production.

The Phoenix 5000 offers 30 second analysis of

Feed, feed ingredients, and by-products

It can test simultaneously for moisture, fat, protein, fiber, starch, ash and more.

If you have already developed databases on other NIR platforms (FOSS[™] Unity Scientific[™]), you can transfer them onto the Phoenix 5000 and begin generating reliable results immediately.

EXPLORE THE PHOENIX

Easy to Use Interface

- 12.3" Microsoft Surface Touchscreen
- Simple tab design
- Fast access to results
- Flexible data handling

All Phoenix instruments are manufactured at Blue Sun Scientific's headquarters in Jessup, Maryland.

NIR Monochromator

The most trusted NIR technology used by the food and agriculture industries



- Scans large and small cups
- Scans entire sample cup
- Integrated standards

Easy to Access Lamp

The 5k-10k hour lamp can be replaced in under 5 minutes with no special tools needed.

Fully Sealed

- Rugged and robust
- Ideal for in-lab or at-line analysis

BLUESUN

HIGHLY ACCURATE, RUGGED AND EASY TO MAINTAIN

Introducing BlueScan

Software Designed for Ease of Use

BlueScan 🛛 🚝 💆 🗙								
	Select Product				Last Sample's Result			
	(SPECTRA ONLY)	CORN	SOYBEAN	BARLEY	Product Name Sample Name			
$\overline{\mathbf{O}}$	WHEAT	SWINE FEED			CONSTITUENT	VALUE	UNIT	
×								
					v. 1.5.4.4 🛛 💋	FLUE SUN		

The Phoenix 5000 is powered by BlueScan software. BlueScan is designed to simplify routine analysis and is packed with features making it ideal for monitoring quality and optimizing production.

BlueScan's tab design ensures your information is just a click away. You can easily view results, historical data, spectra, generate reports and run diagnostics with the press of a button.

Keep Your Database


The Phoenix 5000 Feed Analyzer comes ready to analyze feed ingredients, finished feed and by-products.

If you have existing calibrations that you have built and would like to keep, you can with the Phoenix 5000.

The Phoenix 5000 is backwards compatible with databases from many major NIR manufacturers, including FOSSTM, and Unity ScientificTM.



Monitor and verify incoming ingredients



Control moisture in driers and extruders



Verify finished feed is within specifications

Specifications

Hardware							
Dimensions	14.0 in x 15.0 in x 20.75 in ; 356 mm x 381 mm x 527 mm (LxWxH)						
Weight	34 lbs ; 15.4 kg						
IP Rating	62						
Power Connection	Via external power supply; Input:100-240VAC, 50/60Hz, 65watt Output: 24VDC, 2.7A.						
Light Source	Tungsten Quartz Halogen Lamp, MTBF > 5k-10k hours						
Measurement Mode	Reflectance or Transflectance						
Detector	High Performance InGaAs extended ranged detector						
Wavelength Selection	Scanning diffraction grating monochromator with nominal bandwidth of 10nm (FWHH). Configurable Data presented in 0.5nm, 1nm or 2nm increments						
	Performance						
Wavelength Range	1100-2500 nm						
Number of Data Points	Up to 2800						
Analysis Time	10-60 seconds (24 scans / sample = 32 seconds)						
Wavelength Reproducibility	< 0.02 nm (based on Blue Sun Scientific specific diagnostic test)						
Wavelength Accuracy	< 0.3 nm for Blue Sun selected published peak positions of NIST SRM 1920a standard						
Photometric Noise	<15 µAu in standard range						
User Interface							
Display	12.3" Microsoft Surface Pro Touch Screen running Windows 10						
Data Export Options	LIMS, Modbus, Excel						
	Installation Requirements						
Ambient Temperature	10-40°C						
Storage Temperature	-20-40°C						
Ambient Humidity	<85% RH						
Mechanical Environment	Stationary during use BluePoint						
EMC environment	Laboratory use, Industry requirements TECHNOLOGIES						

Blue Sun Guarantee

Every Phoenix 5000 is backed by our industry leading warranty program which guarantees you receive fast, reliable support from industry experts you can trust at a fair market price. It includes **free**:

- 📞 Phone Support
- Remote Support
- Software Upgrades
- ◯ Calibration Development (1 Year)
- Warranty on Parts (1 Year)